

FloraQUICK®



Range of application:

FloraQUICK® is an alcohol-based ready to use spray or wipe-on disinfectant for quick disinfecting and cleaning of medical products, surfaces, and objects found in a clinic.

The extensive scope of protection and fast-acting effect enables a versatile usage of FloraQUICK® for purposes of safe sterilization.

Expertises:

Prof. Dr. med. H.-P. Werner, A-Bischofshofen, DGHM-procedure, 03.2005

Prof. Dr. med. H.-P. Werner, A-Bischofshofen, TbB-effi acy, 03.2005

Dr. med. habil. Schrader, Weimar, DGHM-procedure, 04.2005

Dr. J. Steinmann, Bremen, virucidal effi acy, 05.2005

Prof. Dr. med. H.-P. Werner, A-Bischofshofen, bactericidal effi acy (strains of MRSA), 03.2005

Ready to use quick disinfectant for medical devices and surfaces

- Fully disinfects in 30 seconds
- No phenols or aldehydes
- Pleasant smell
- Ready to use

Extensive efficacy spectrum:

- Bactericidal incl. strains of MRSA
- Tuberculocidal
- Levurocidal
- Limited virucidal (inactivates enveloped viruses such as HBV/HCV/HIV)

Active ingredients:

100 g of FloraQuick® contain:

37.0 % 1-propanol

24.0 % ethanol

Exposure time:

Bacteria (incl. TbB and strains of MRSA) as per DGHM¹/VAH², fungus (*Candida albicans*) and viruses (incl. HBV, HCV, HIV):

High organic load 1 min.

Low organic load 30 sec.

Unit	PZN	Article-No.
250 ml aerosol can	0260741	111
1 L metered bottle	0260787	112
Sprayhead for 1L bottle		003
5 L canister	0281134	113

On request we would be pleased to submit tailored conditions.

+49-721-62686-0

info@favodent.com

¹ German Society for Hygiene and Microbiology

² Association for Applied Hygiene