

Floracin® NF Concentrate



Range of application:

Floracin® NF is an alcoholfree concentrate that both disinfects and cleans large washable surfaces (walls and floors), surfaces of objects and medical products in clinics, doctors' practices and laboratories. It is also suitable for disinfecting prosthetics and silicon impressions. Due to its excellent grease dissolving properties Floracin® NF Concentrate can also be used in the food sector.

This quaternary bondbased product has a broad range of effectiveness and is fastacting and does not cause any damage to materials - even with high concentrations of germs.

Disinfectant for cleaning and disinfection for surfaces of medical devices

- No alcohols or aldehydes
- Gentle and compatible with materials
- Very good cleaning properties
- Excellent grease dissolving properties
- Biodegradable

Extensive efficacy spectrum:

- Bactericidal
- TbB (Mycobact. terrae und M. avium)
- Limited virucidal (inactivates enveloped viruses such as HBV/ HIV/ HCV)
- Levurocidal

Active ingredients:

100 g Floracin® NF Concentrate contain:
3,15g Didecylmethylpolyoxyethylaminopropionate
3,5g Bis (3-aminopropyl)dodecylamine, surfactants

Exposure time:

Bakteria und Candida albicans acc. DGHM/VAH (high organic load)	1,0 % - 5 min 0,5 % - 15 min
Enveloped viruses (incl. HBV, HIV, HCV) acc. BGA(RKI)/DVV	1,0 % - 5 min 0,5 % - 15 min
TbB (Mycobact. terrae) (without organic load)	1,0 % - 30 min 0,5 % - 60 min
TbB (M. terrae und M. avium)	1,5 % - 60 min

Unit	Article-No.
1 l metered bottle	155
5 l canister	156

On request we would be pleased to submit tailored conditions.

+49 (0) 721/62 6860
info@favodent.com

Surface Disinfection

N-56326

€ 0482