

## Temporary sealing compound

- Ready to use
- Easy-to-mold consistency
- Excellent adhesion
- Strong border sealing
- Eugenol-free
- Cure time of only 30 min

### Composition:

100 g of Provis® contains:

37.2 g zinc oxide  
10.4 g potassium sulfate  
9.6 g zinc sulfate • 1H<sub>2</sub>O

**Color:** white

Unit	Article-No.
40 g Jar	601
3 x 40 g Economy Pack	602

On request we would be pleased to submit tailored conditions.

+49-721-62686-0  
info@favodent.com

CE 0482

# Provis®



### Range of application:

Provis® is a temporary zinc oxide-based sealing compound for cavities. It has a particularly smooth and pliable consistency, optimal border sealing and a low cure time only 30 minutes.

Adhesion is excellent even on moist dentin. This product is always easy to mold.

Provis® does not irritate the dental pulp and can be easily removed with a probe or an excavator.



Dental Supplies

REF 6018-02 2018-09 EN