

## Versatile device for dentists

- Pharmaceutical quality
- Indispensable for each practice and clinic
- Multi-purpose

### Composition:

Vaselinum album - DAB product

Color: pure white

Unit	Art.-Nr.
20 g Tube	720
12 x 20 g Tube	712
30 g Can	730
200 g Can	700

On request we would be pleased to submit tailored conditions.

+49-721-62686-0  
info@favodent.com



# Vaseline



### Range of application:

Vaseline has many applications in dentistry, for example:

- Preparing the buccal area prior to taking impressions
- Mixing with cements based on zinc oxide eugenol when trying out crowns and bridges to ease their removal
- Greasing a bridge member (protecting the crown opening) to facilitate the removal of residual cement after hardening
- Covering composite filling and/or veneers before putting acrylate on bite correction braces and similar appliances to prevent polymerization
- Isolating matrices when dealing with glass ionomer filling



Dental Supplies

REF 7208-01 2016-06 EN